



Construct, Conduct & Comprehend Physics Experiments

## LEEI-62 Interference, Diffraction & Polarization of Microwave - Advanced

- Serial experiments
- Quantitative measurement
- Easy to use, stable and reliable
- Affordable





Bragg diffraction



Microwave and lightwave are both electromagnetic waves. They share common phenomena of all waves, such as reflection, refraction, polarization, interference, and diffraction. However, as the wavelength of a microwave is about 4 orders of magnitude larger than that of a visible lightwave, the experimental phenomena and apparatus of microwave are different.

This LEEI-62 experimental system of microwave interference, diffraction and polarization consists of a microwave transmitter with a variable attenuator, a microwave detector with amplifier and meter display, a transmitting horn and a receiving horn, as well as other accessories such as single slit, double slits, beam splitter, crystal model, etc. It can provide quantitative measurements of a series of experiments as:

- 1. Verify reflection law of electromagnetic wave
- 2. Measure intensity distribution of single-slit diffraction
- 3. Measure intensity distribution of double-slit interference
- 4. Measure microwave wavelength through interference of Michelson interferometer
- 5. Study polarization and verify Malus' law
- 6. Verify Bragg diffraction of crystal model





## Specifications

Microwave frequency	Range: 8.6~9.6 GHz; draft: ± 5×10 <sup>-4</sup> /15-min; error: ±40 MHz
Microwave output power	>20 mW
Working voltage	DC 12 V
Wave form	Equal amplitude
Internal modulation	Square-wave at 1 kHz
Output waveguide dimensions	Inner: 22.86 mm × 10.16 mm
Standing-wave coefficient	≤ 1.2 (synthesis voltage)
Output microwave wavelength	λ =32.02 mm (factory preset)
Horn antenna	Gain ~ 20 dB; lobe width H 20 <sup>0</sup> & E 16 <sup>0</sup>

## **Part List**

3-cm solid state microwave signal source	
Base	
Indexing disk	1
Fixed arm	1
Horn antenna	2
Variable attenuator	1
Crystal detector	1
Indicator meter	1
Video cable	1
Reflection plate	1
Single-slit plate (slit width=70 mm)	
Double-slit plate (slit width=40 mm; separation=170 mm)	
Partial transparent plate	1
Simulated crystal with support	1
Reading mechanism	1
Support base	1
Support post	
Module piece	1
Screws and nuts (set)	1





Lambda Scientific Systems, Inc. 16300 SW 137th Ave, Unit 132

Miami, FL 33177, USA

Phone: 305.252.3838 Fax: 305.517.3739

E-mail: sales@lambdasys.com Web: www.lambdasys.com

Note: above product information is subject to change without notice.