

Construct, Conduct & Comprehend Physics Experiments

LLD-1 VGA Video Camera with Length Measurement Function



- No computer and software needed
- 1.3 megapixel, color, high quality CMOS, resolution 1280 x 1024, clarity > 650 TV lines
- Up to 8 crosshair markers movable to any location on screen
- Length measurement in both H & V directions after calibration
- High measurement accuracy up to sub-micron
- Live VGA video output to display, projector or digital TV
- Easy-to-use by menu operation using keypad with clear & stable image

Applications

- Academic teaching
- Microscopic amplification systems
- QC/QA in micro IC assembly
- PCB and SMT production
- Industrial precision positioning
- Medical and scientific R&D
- Jewelry appraisal

Installation

- Install a proper camera lens onto the CS mount of the camera
- Connect the VGA OUT terminal to a VGA display using a VGA cable
- Plug the 5V DC AC adaptor into the 5V DC IN socket of the camera

Menu Operation

- Press MENU once, main menu pops up; press MENU again, main menu closes.
- Press \uparrow/\downarrow to highlight menu, then press \rightarrow to select menu.
- If there is a submenu, press ↑/↓ to highlight submenu, then press ←/→ to change value of the selected submenu.
- Select Exit to return to main menu.



Method of Measurement


- Place the sample under test on the specimen stage, while keeping the calibration of the system unchanged.
- Through keys SEL,↑,↓,← and → to move lines to measure distance between a pair of lines in H or V direction
- Screen freeze: press POWER for one second, screen will freeze for more stable measurements.
- Press POWER one second again to release freeze function.
- Press POWER about 3 seconds to turn off power of the camera.

Specifications

Sensor	1.3 megapixel, color, ¼" CMOS
Effective pixels	1280 (H) x 1024 (V), pixel size 2.8 μm x 2.8 μm
Data	8 bits
Frame rate	15 frames/s progressive scan, electronic shutter
Clarity	> 650 TV lines
Dynamic range	> 63 dB
S/N	44 dB
Sensitivity	1 V/lux-sec (at 550 nm)
Exposure	10 ms ~ 60 ms, auto/manual
White balance	auto
Lens mount	CS/C
Video output	VGA, resolution selectable1024 x 768, 1280 x 1024
Menu function	brightness, contrast, sharpness, saturation, blue & red offset, auto exposure
	H-mirror, V-mirror, negative, crosshair number, length calibration
Working temperature	0~60 °C
Power supply	5 VDC ± 10%, 400 mA
Dimensions	50 x 50 x 90 (mm)
Package	camera unit (CS), 5 VDC adapter, keypad, manual



Microscopic PCB inspection



Newton's ring experiment



Michelson interferometer

Lambda Scientific Systems, Inc. 16300 SW 137th Ave, Unit 132 Miami, FL 33177, USA

Phone: 305.252.3838 Fax: 305.517.3739

E-mail: sales@lambdasys.com Web: www.lambdasys.com

Note: above product information is subject to change without notice.