

3.5 Constructing a Mach-Zehnder Interferometer

Objective:

Learn to construct a Mach-Zehnder interferometer and observe the interference phenomena

Experimental Setup



Figure 5-1 Photo of experimental setup




Figure 5-2 Configuration of components

Note: Photos may vary from actual parts

- | | |
|---|----------------------------------|
| 1: He-Ne Laser L (LLL-2) | 11: Two-axis Tilt Holder (SZ-07) |
| 2: Laser Holder (SZ-42) | 12: Carrier (LEPO-54-2) |
| 3: X-Trans Carrier (LEPO-54-3) | 13: Magnetic Base (SZ-04) |
| 4: Beam Expander ($f'=15$ mm) | 14: White Screen (SZ-13) |
| 5: Lens Holder (SZ-08) | 15: Beam Splitter (5:5) |
| 6: X-Z-trans Carrier (LEPO-54-4) | 16: Lens Holder (SZ-08) |
| 7: Beam Splitter (5:5) | 17: Magnetic Base (SZ-04) |
| 8: Two-axis Tilt Holder (SZ-07) or Plate Holder (SZ-12) | 18: Flat Mirror |
| 9: Carrier (LEPO-54-2) | 19: Two-axis Tilt Holder (SZ-07) |
| 10: Flat Mirror | 20: Magnetic Base (SZ-04) |

Principle

Mach-Zehnder interferometer (MZI) is also derived from Michelson interferometer. It is used for studying the change in the wave front when the light wave passing through an object of interest. The schematic diagram of the Mach-Zehnder interferometer is shown in Figure 5-3. The light beam from a He-Ne laser is first divided by a beam splitter BS_1 into two beams of equal intensity. After reflected by mirrors M_1 and M_2 , the two beams are recombined via a second beam splitter BS_2 . Then the interference pattern will be observed on a view screen.

A MZI is frequently used in the fields of plasma physics, aerodynamics, and heat transfer for the measurement of density, pressure, and temperature changes in gases.




Figure 5-3 Schematic of Mach-Zehnder interferometer

Experimental Procedure

1. Refer to Figure 5-2, align all components in same height on the optical rail. At the moment, beam expander should not be included in the light path;
2. Adjust the output of the He-Ne laser to make it parallel to experimental table (and along the optical rail if applicable);
3. Adjust beam splitter BS_1 at an angle of 45° with respect to the beam axis, and adjust its tilt so that the two beams (transmission and reflection) are parallel to the table;
4. Adjust mirrors M_1 and M_2 until the light beams reflected by them strike at same position on BS_2 ;
5. Put in the beam expander, interference pattern will be observed on the screen;
6. Clap the experimental table and see any changes of the interference pattern.