

LEOK-3-12 Studying on Young's Double-Slit Interference

- Complete set
- Cost effective solution
- Detailed instructional manual
- Easy alignment





Figure 12-1 Schematic of experiment setup

- | | |
|--------------------------------------|------------------------------|
| 1: Sodium Lamp (LLE-2, w/small hole) | 7: Biprism Holder (SZ-41) |
| 2: Lens L_1 ($f' = 50$ mm) | 8: Double-slit Plate D |
| 3: Lens Holder (SZ-08) | 9: Adapter Piece (SZ-09) |
| 4: Adjustable Slit S (SZ-27) | 10: DMM Holder (SZ-36) |
| 5: Lens L_2 ($f' = 150$ mm) | 11: Eyepiece of DMM |
| 6: Lens Holder (SZ-08) | 12-16: Magnetic Base (SZ-04) |

Theory

To get an interference pattern, the two beams leaving from the slits must have same frequency and a definite phase relation. Generally speaking, most light sources cannot satisfy this condition. In 1801, Thomas Young allowed a single, narrow beam of light to fall on two narrow, closely spaced slits. He placed a viewing screen opposite to the slits. When the light from the two slits struck the screen, a regular pattern with alternative dark and bright rings appeared. When first performed, Young's experiment offered an important evidence for the wave nature of light. The schematic of Young's double-slit interference is shown in Figure 12-2.

In this way, the light emitted from S_2 and S_3 has a definite phase relation because the secondary wave sources from the same wave surface S_1 are always coherent. The light path difference (d is the distance between the two slits of the double-slit plate) is:



$$\delta = r_2 - r_1 \approx d \sin \theta \approx d \tan \theta = d \frac{x}{D} \quad (12-1)$$

Where D is the distance between the viewing screen and the slits, x is the vertical distance between the viewing location and the center of the double slits, and θ is a half of the viewing angle between the lines from the viewing point on the screen to the two slits.

Figure 12-2 Schematic of Young's double-slit experiment

If the path difference between a particular point on the screen to the two slits is equivalent to a half of the wavelength (or multiples thereof) of the light, then complete destructive interference will occur at that point, and thus a dark spot will be observed.

$$\delta = d \frac{x}{D} = \pm(2k + 1) \frac{\lambda}{2} \quad \text{(Dark interference fringes)} \quad (12-2)$$

Conversely, if the path difference is equivalent to an integer multiple of the wavelength of the light, then complete constructive interference will occur, and a bright spot will appear on the screen.

$$\delta = d \frac{x}{D} = \pm k \lambda \quad \text{(Bright interference fringes)} \quad (12-3)$$

So the distance between two adjacent dark fringes (or bright fringes) is:

$$\Delta x = \frac{D}{d} \lambda \quad (12-4)$$

In this formula Δx and D can be measured, so if we know one of d and λ , we can calculate the other. In this experiment, if a laser rather than a Sodium lamp is used as the source, the experiment will be easier to conduct and the interference fringes will be observed more obviously.

Experiment Procedures

1. Refer to Figure 12-1, align all components in same height along a straight line;
2. Focus the small hole of the light source onto the single slit S by lens L_1 (i.e. collect as more as possible light power to the single slit);
3. Remove the double-slit D temporary, adjust S to relatively wide (e.g. 0.5 mm), use lens L_2 to image single slit S onto the scale plane of the eyepiece M (note: the distance between S and M must be larger than 4 times focal length of L_2).
4. Place back the double-slit D and let it close to L_2 , look into the eyepiece M, narrow single slit S, finely rotate D to make S parallel to D in vertical direction strictly. Double-slit interference pattern with equal-interval bright/dark pairs will be observed;
5. Further finely adjust single slit width and distances between components to optimize interference fringe.
6. Measure the interval e between two adjacent fringes using direct measurement microscope, also measure the distance L between double-slit plate and the microscope;
7. Use the known value of double-slit interval t and expression of

$$e = \frac{L \lambda}{t}, \text{ so that the wavelength } \lambda \text{ of the illumination light can be obtained.}$$

