

## LEOK-3-9 Build a Microscope and Determine Magnification Power

- Complete set
- Cost effective solution
- Detailed instructional manual
- Easy alignment




Figure 9-1 Schematic of experiment setup

- |                                      |   |
|--------------------------------------|---|
| 1: Sodium Lamp S (LLE-2)             | 9: Millimetre Ruler $M_2$ ( $l = 30$ mm)                                |
| 2: Reticle $M_1$ (1/10 mm)           | 10: Bi-prism Holder (SZ-41)   |
| 3: Kinematic Holder (SZ-07)          | 11, 14: Z-adjustable Stage (SZ-03)                                      |
| 4: Objective Lens $L_o$ ( $f=45$ mm) | 12: Three-Axis Stage (SZ-01)  |
| 5,6: Lens Holder (SZ-08)             | 13: Two-Axis Stage (SZ-02)  |
| 7: Eyepiece Lens $L_e$ ( $f=29$ mm)  | Other: White light source (LLC-3) for illuminating the millimetre ruler |
| 8: 45° Glass Holder (SZ-45)          |   |

### Theory

As shown in Figure 9-2, the optical system of a microscope employs an objective with a short focal length and a magnifying eyepiece. The magnification is achieved in two stages as shown in Figure 9-2. The microscope objective forms an enlarged image of the object in a position that is suitable for viewing through the eyepiece; the magnification through the objective is given by

$$y_2/y_1 = \Delta / f_o'$$

Generally speaking, the focal length of the eyepiece  $f_e'$  is much less than the distance of the image from the eyepiece  $D$ , (for normal sight,  $D$  is approximate 250 mm), so

$$y_3/y_2 \approx D / f_e'$$

Then we get the total magnification:

$$M = \frac{y_3}{y_1} = \frac{y_3}{y_2} \frac{y_2}{y_1} = \frac{D\Delta}{f_o' f_e'}$$

Where  $\Delta$  is the distance between the focus of objective and the focus of eyepiece,  $f_o'$  is the focal length of objective and  $f_e'$  is that of eyepiece.




Figure 9-2 Schematic of microscope imaging

## Experiment Procedures

1. Refer to Figures 9-1, align all components in same height;
2. Fix the spacing between  $L_o$  and  $L_e$  at 240 mm;
3. Move reticle plate  $M_1$  back and forth, till a clear  $M_1$  virtual image is observed behind  $L_e$ ;
4. Mount the 45° Glass Holder onto the post of the lens holder of the eyepiece lens  $L_e$ , set the glass surface at 45° angle with respect to the optical axis; the glass is acting as a beam splitter (BS);
5. Put the millimetre ruler  $M_2$  beside the Glass Holder (perpendicular to the main optical axis of the microscope system) and approximate 250 mm distance from the 45° glass; place the white light source (LLC-3) behind of the millimetre ruler to illuminate the ruler with proper intensity;
6. View behind the Glass by one eye, finely rotate the Glass holder to overlap the microscope virtual image from  $M_1$  and the  $M_2$  image from the glass reflection; note: if the brightness difference between the two images is too much, try to adjust the intensity of LLC-3;
7. Finely adjust  $M_1$  to eliminate viewing difference between the two images;
8. Count the scale amount  $a$  in  $M_1$  image included in the range of 30 mm of image  $M_2$ ;
9. Calculate the measured magnification of the assembled microscope and its theoretical magnification:

Measured Magnification: 
$$M = \frac{30 \times 10}{a}$$

Theoretical Magnification: 
$$M' = \frac{25\Delta}{f_o' f_e'}$$

where, 
$$\Delta = D - (f_o' + f_e')$$

