

4.4 Constructing a Slide Projector

Components Required: white light source S, condensing lens L_1 ($f=50$ mm), lens holder, transparent slide P, plate holder, projection lens L_2 ($f=190$ mm), 2-axis tiltable holder, white screen.

Principle

As shown in Figure 4.4-1, L_1 is a condenser, and L_2 is a projection lens. A slide is just behind L_1 (we can assume $v_1 = u_2$). If the magnification of a slide projector is M , the length of slide projector is D , and the focal length of L_1 and L_2 are f_1 and f_2 , respectively.

By taking $M = v_2/u_2$, $1/f_2 = 1/u_2 + 1/v_2$, we can get

$$f_2 = \frac{1}{M+1} v_2 \quad (4.4-1)$$

By taking $D = u_1 + v_1$, $v_1 = u_2$, $1/f_1 = 1/u_1 + 1/v_1$, we can get

$$f_1 = \frac{v_2}{M} - \frac{1}{D} \left(\frac{v_2}{M} \right)^2 \quad (4.4-2)$$




Figure 4.4-1 Schematic of slide projector imaging




Figure 4.4-2 Configuration of components

Experimental Procedures:

- 1) Refer to Figure 4.4-2, align all components in same height on the optical rail, set the distance between L_2 and screen H about 0.8 m (if space available, it should be better to set a larger distance);
- 2) Move slide P back and forth, till a clear image (imaged by L_2) is observed on H ;
- 3) Fix condenser L_1 as close as possible to P (may use the adapter piece SZ-09), remove P , move light source S back and forth, till the image of S formed by L_1 is clear on L_2 aperture plane;
- 4) Put back slide P at its pervious location, observe the brightness and uniformity of the projected image on the screen;
- 5) Remove L_1 , observe the brightness and uniformity of the projected image again, and recognize the function of L_1 .
- 6) Through the experiment, understand the principle of a slide projector and the function of its condenser, learn how to adjust a projection optical system, and understand illuminating condition for achieving a uniform light field on the screen (Kohler illumination).