

4.3 Determining the magnifications of a microscope and a telescope

Components required: small illumination lamp, white light source, plate holder, 1/10 mm reticle, millimetre ruler, biprism holder, lens holder(2), 45° glass holder, microscope eyepiece, stand ruler, adapter piece, and lenses $f=45$ and 190 mm.

4.3.1 Magnification of a microscope

As shown in Figure 4.3-1, the optical system of a microscope employs an objective with a short focal length and a magnifying eyepiece. The magnification is achieved in two stages as shown in Figure 4.3-1. The microscope objective forms an enlarged image of the object in a position that is suitable for viewing through the eyepiece; the magnification through the objective is given by

$$y_2/y_1 = \Delta / f_o' \quad (4.3-1)$$

Generally speaking, the focal length of the eyepiece f_e' is much less than the distance of the image from the eyepiece D , (for normal sight, D is approximate 250 mm), so

$$y_3/y_2 \approx D/f_e' \quad (4.3-2)$$

Then we get the total magnification:

$$M = \frac{y_3}{y_1} = \frac{y_3}{y_2} \frac{y_2}{y_1} = \frac{D\Delta}{f_o'f_e'} \quad (4.3-3)$$

Where Δ is the distance between the focus of objective and the focus of eyepiece, f_o' is the focal length of objective and f_e' is that of eyepiece.




Figure 4.3-1 Schematic of microscope imaging




Figure 4.3-2 Schematic of experiment setup

- | | |
|---|---|
| 1: Bromine Tungsten Lamp S (LLC-3) | 8: 45° Glass Holder (SZ-45) |
| 2: Reticule M_1 (1/10 mm) | 9: Optical Rail (LEPO-54) |
| 3: Biprism Holder (SZ-41) | 10: Magnetic Base (SZ-04) |
| 4: Objective Lens L_o ($f'_o=34$ mm) | 11: Lens Holder (SZ-08) |
| 5: Two-axis Tilt Holder (SZ-07) | 12: Millimetre Ruler M_2 ($l=30$ mm) |
| 6: Eyepiece Lens L_e ($f'_e=45$ mm) | Other: Small illuminating lamp (LLC-6) mounted on plate holder SZ-12, adapter piece SZ-09A. |
| 7: Lens Holder (SZ-08) | |

Experiment Procedures:

- 1) Refer to Figures 4.3-2, align all components in same height along the rail;
- 2) Fix interval between L_o and L_e as $D = 240$ mm;
- 3) Move reticle plate M_1 back and forth, till a clear M_1 virtual image is observed behind L_e ;
- 4) Mount the beam splitter (B) behind the lens L_e on the post of lens holder 8 (the glass plate is at 45° angle with respect to the optical axis);
- 5) Place the millimetre ruler M_2 beside the B (vertical to main optical axis along the rail) and approximate 250 mm distance from B (*in diagram*); it is illuminated by the small illuminating lamp from its back; (i.e. mount the millimetre ruler M_2 on a lens holder; then install the millimetre ruler and the small illuminating lamp at the two ends of the adapter piece SZ-09A.)
- 6) View behind B by one eye, finely rotate B 's angle to overlap the microscope virtual image from M_1 and the M_2 image from the glass reflection; adjust the brightness of the LLC-3 lamp to achieve close brightness of the two images.
- 7) Finely adjust M_1 to eliminate viewing difference between the two images;
- 8) Count the scale amount a in M_1 image included in the range of 30 mm of image M_2 ;

- 9) Calculate the measured magnification of the assembled microscope and its theoretical magnification:

$$\text{Measured Magnification: } M = \frac{30 \times 10}{a}$$

$$\text{Theoretical Magnification: } M' = \frac{250 \times \Delta}{f_o' f_e'}, \text{ where, } \Delta = D - (f_o' + f_e').$$

4.3.2 Magnification of a telescope

As seen in Figure 4.3-3, the magnifying power of a telescope used for observing an object at infinity is defined as the angular magnification at the pupil because the angles are very small:

$$M = \frac{\tan \omega'}{\tan \omega} = \frac{\omega'}{\omega} = \frac{f_o'}{f_e'} \quad (4.3-4)$$

where f_o' and f_e' are the focal lengths of the objective and eyepiece lenses, respectively, and ω' and ω are the object and image angles at the eyepiece and objective lenses, respectively.




Figure 4.3-3 Schematic of telescope imaging at infinity

As shown in Figure 4.3-4, when observing an object at quasi-infinity, the power of magnification is:

$$M = \frac{\tan \omega'}{\tan \omega} = \frac{y_2 / s_2}{y_1 / (s_1 + s_1' + s_2)} \quad (4.3-5)$$

Since $y_2 / y_1 = s_1' / s_1$, therefore,

$$M = s_1' (s_1 + s_1' + s_2) / s_1 s_2 \quad (4.3-6)$$




Figure 4.3-4 Schematic of telescope imaging at quasi-infinity.

Experiment Procedures:

- 1) Refer to Figure 4.3-4, align L_o and L_e in same height on the optical rail with their spacing about f_o+f_e , place the stand ruler in front of L_o at a distance of about 3 meters;
- 2) Move objective lens L_o back and forth, behind L_e , use one eye to observe the image of the ruler till it is clear;
- 3) Use another eye to directly observe the scale marks on the ruler, count how many scale marks (amount M') on the ruler image directly to the eye are contained in one mark of the magnified image of the telescope; M' is the measured magnification of the telescope;
- 4) Calculate magnification M using Eq. (4.3-4), compare it with M' .



4.4 Constructing a Slide Projector

Components Required: white light source S, condensing lens L_1 ($f=50$ mm), lens holder, transparent slide P, plate holder, projection lens L_2 ($f=190$ mm), 2-axis tiltable holder, white screen.

Principle

As shown in Figure 4.4-1, L_1 is a condenser, and L_2 is a projection lens. A slide is just behind L_1 (we can assume $v_1 = u_2$). If the magnification of a slide projector is M , the length of slide projector is D , and the focal length of L_1 and L_2 are f_1 and f_2 , respectively.

By taking $M = v_2/u_2$, $1/f_2 = 1/u_2 + 1/v_2$, we can get

$$f_2 = \frac{1}{M+1} v_2 \quad (4.4-1)$$

By taking $D = u_1 + v_1$, $v_1 = u_2$, $1/f_1 = 1/u_1 + 1/v_1$, we can get

$$f_1 = \frac{v_2}{M} - \frac{1}{D} \left(\frac{v_2}{M}\right)^2 \quad (4.4-2)$$




Figure 4.4-1 Schematic of slide projector imaging




Figure 4.4-2 Configuration of components

Experimental Procedures:

- 1) Refer to Figure 4.4-2, align all components in same height on the optical rail, set the distance between L_2 and screen H about 0.8 m (if space available, it should be better to set a larger distance);
- 2) Move slide P back and forth, till a clear image (imaged by L_2) is observed on H ;
- 3) Fix condenser L_1 as close as possible to P (may use the adaptor piece SZ-09), remove P , move light source S back and forth, till the image of S formed by L_1 is clear on L_2 aperture plane;
- 4) Put back slide P at its pervious location, observe the brightness and uniformity of the projected image on the screen;

- 5) Remove L_1 , observe the brightness and uniformity of the projected image again, and recognize the function of L_1 .
- 6) Through the experiment, understand the principle of a slide projector and the function of its condenser, learn how to adjust a projection optical system, and understand illuminating condition for achieving a uniform light field on the screen (Kohler illumination).

